Gynecare **PROSIMA**TM

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Manufactured for:
ETHICON Women's Health & Urology
A division of ETHICON, INC.
a Johnson & Johnson company
Somerville, New Jersey 08876-0151
Made in Switzerland
© Ethicon, Inc. 2008

Legal Manufacturer
ETHICON, S.r.l.
Rue du Puits-Godet 20
CH-2000 Neuchâtel
Switzerland
P21070/D
STATUS: 02/2010
LAB0011038.2

EXHIBIT E

ETH.MESH.02341398


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7


Figure 8A


Figure 8B


Figure 8C


Figure 8D


Figure 8E


Figure 8F


Figure 9A


Figure 9B


Figure 9C


Figure 9D


Figure 9E


Figure 10


Figure 11


Figure 12


Figure 13


Figure 14


Figure 15


Figure 16


Figure 17


Figure 18


Figure 19

ENGLISH

A – Mesh Implant in Implant Carrier
 B – Vaginal Support Device (VSD)
 C – Balloon
 D – Anterior Insertor
 E – Posterior Insertor
 F – Syringe
 G – Straps of Mesh Implant
 H – Pockets of Mesh Implant
 I – Apical tab of Mesh Implant
 J – Distal groove of Mesh Implant
 K – Body of Mesh Implant
 L – Connection groove of Posterior Insertor
 M – Distal end of Posterior Insertor
 N – Trimable sections of VSD
 O – Suture eyelets of VSD
 P – Balloon
 Q – Connector plug of Balloon
 R – Valve of Balloon
 S – Cap of Balloon
 T – Inflation line of Balloon

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ENGLISH

Please read all information carefully.

Failure to properly follow instructions may result in improper functioning of the devices and lead to injury.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Training on the use of the GYNECARE PROSIMA™ Pelvic Floor Repair Systems is recommended and available. Contact your company sales representative to arrange for this training.

INDICATIONS

The GYNECARE PROSIMA Pelvic Floor Repair Systems, through the placement of GYNECARE GYNEMESH™ PS Nonabsorbable PROLENE™ Soft Mesh Implants, are indicated for tissue reinforcement and long-lasting stabilization of fascial structures of the pelvic floor, either as mechanical support or bridging material for the fascial defect. The Systems provide maintenance of the vaginal canal during the period of healing following surgical repair of vaginal wall prolapse, while supporting the position of the Mesh Implants.

DESCRIPTION

The GYNECARE PROSIMA Anterior, Posterior, and Combined Pelvic Floor Repair Systems consist of pre-cut GYNECARE GYNEMESH PS Mesh Implant(s), and instruments to facilitate Mesh Implant placement and postoperative support (see Figure 1). The following table summarizes the components included with each System:

PELVIC FLOOR REPAIR SYSTEM	COMPONENTS (see Figure 1)				
	Mesh Implant in Carrier (A)	Vaginal Support Device – Balloon Assembly (B&C)	Anterior Insertor (D)	Posterior Insertor (E)	Syringe (F)
Anterior	1	1	1		1
Posterior	1	1		1	1
Combined	2	1	1	1	1

Table 1 – GYNECARE PROSIMA Pelvic Floor Repair System Components

GYNECARE GYNEMESH PS

GYNECARE GYNEMESH PS is mesh constructed of knitted filaments of extruded polypropylene identical in composition to PROLENE™ Polypropylene Suture, (ETHICON, INC.). This material, when used as a suture, has been reported to be nonreactive and to retain its strength indefinitely in clinical use. The mesh affords excellent strength, durability, and surgical adaptability, with sufficient porosity for necessary tissue ingrowth. Blue PROLENE Suture monofilaments have been incorporated to produce contrast striping in the mesh. The mesh is constructed of reduced-diameter monofilament fibers, knitted into a unique design that results in a mesh that is approximately 50 percent more flexible than standard PROLENE™ Polypropylene Mesh. The mesh is knitted by a process which interlinks each fiber junction and which provides for elasticity in both directions. This construction permits the mesh to be cut into any desired shape or size without unraveling. The bi-directional elastic property allows adaptation to various stresses encountered in the body.

Mesh Implant

The Mesh Implant is constructed from GYNECARE GYNEMESH PS. The Mesh Implants are pre-cut in a Y-shape for repair of anterior, posterior and/or apical vaginal defects. See Figure 2. The Mesh Implant has 2 straps and a central body. There is an apical tab on the proximal end for tacking down with suture to minimize movement of the Mesh Implant during strap placement. There is a distal groove on the distal end to aid in alignment of the Mesh Implant. There are pre-formed pockets on the Mesh Implant straps to enable placement with the Inserters. The Mesh Implant is provided in an Implant Carrier comprised of uncoated Tyvek® and a clear plastic film, designed for easy removal of the Mesh Implant.

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The Anterior Insertor is a single-patient-use instrument designed to facilitate the insertion of the Mesh Implant straps into the previously dissected anterior tissue channels. **NOTE: The Anterior Insertor is not intended to dissect tissue.** The Anterior Insertor is designed to be compatible with the Mesh Implant pockets to enable strap placement on both sides of the patient in the anterior compartment. See Figures 3 and 4.

Posterior Insertor

The Posterior Insertor is a single-patient-use instrument designed to facilitate insertion of the Mesh Implant straps into the previously dissected posterior tissue channels. **NOTE: The Posterior Insertor is not intended to dissect tissue.** A standard needle holder/driver attaches to the Posterior Insertor as a stabilizer for controlled insertion. The Posterior Insertor is designed to be compatible with the Mesh Implant pockets to enable strap placement on both sides of the patient in the posterior compartment. See Figure 5.

Vaginal Support Device (VSD)

The VSD is a single-patient-use device designed to provide postoperative support for the vaginal tissues after placement of the mesh and closure of the vaginal incision(s). The apical end is the widest end of the VSD and contains trimmable sections. After initial sizing in the patient, the size of the VSD can be adjusted to conform to the patient's anatomy by trimming off designated apical sections. The VSD resides in the upper 2/3 of the vagina for 3 to 4 weeks and is then removed from the patient. See Figure 6.

Balloon

The Balloon is a single-patient-use device designed to replace postsurgical vaginal gauze packing. The Balloon's volume is adjustable in order to fill the vaginal canal, and for abutment of the vaginal wall to the Mesh Implant. The Balloon is provided pre-attached onto the VSD. Figure 7 shows the deflated Balloon without the VSD attached. The Balloon remains in the patient for up to 1 day.

Syringe

A 50-mL syringe is provided to inflate the Balloon.

SECTION 1: PRINCIPLES OF THE PROCEDURE USING THE GYNECARE PROSIMA SYSTEM

A pelvic floor repair procedure using the GYNECARE PROSIMA System aims to achieve an anatomical, durable, and standardized repair of pelvic organ prolapse. Depending on the site of the prolapse and surgeon's preference, the repair can be anterior and/or posterior. Hysterectomy or uterine conservation can be combined with the procedure using the GYNECARE PROSIMA System. If indicated, a perineal repair or a suburethral sling for the treatment of stress urinary incontinence can be performed concomitantly when using the GYNECARE PROSIMA System. A retropubic or transobturator suburethral sling can be used.

The prolapse repair is achieved by the placement of 1 or 2 Mesh Implants via a vaginal approach. At the completion of surgery, a VSD with an inflatable Balloon is placed in the vagina for sizing and then the VSD is sutured into place, thus supporting the vagina and the Mesh Implant(s) during tissue ingrowth. Once inflated, the Balloon replaces traditional gauze packing by filling the vaginal cavity and abutting the Mesh Implant(s) to the vagina. The day after surgery, the Balloon is deflated and removed from the vagina without dislodging the VSD. The VSD remains in place for a maximum of 4 weeks following surgery, during tissue ingrowth into the Mesh Implant(s).

SECTION 2: RATIONALE FOR THE GYNECARE PROSIMA SYSTEM

Following conventional surgery for pelvic organ prolapse, the repaired tissues are exposed to increases in intra-abdominal pressure as the patient mobilizes, coughs, vomits, and strains with bowel evacuation. These rises in intra-abdominal pressure may adversely affect the healing of the vaginal repair and may lead to surgical failure and recurrent prolapse. By reinforcing the vaginal repair with the Mesh Implant and supporting the vagina with the VSD for 3 to 4 weeks following surgery, the GYNECARE PROSIMA System is designed to reduce the risk of surgical failure and recurrent prolapse.

During the anterior vaginal repair the body of the Mesh Implant is intended to be placed without tension between the urinary bladder and the upper 2/3 of the vagina, extending laterally at the level of the arcus tendineus fascia pelvis (ATFP). During the posterior vaginal repair, the Mesh Implant body is intended to be placed without tension between the rectum and upper 2/3 of the vagina, fitting laterally over the levator ani muscles. The apical section of the Mesh Implant body is intended to reach the vaginal apex. Anteriorly, the Mesh Implant may be tacked to the pre-vesical tissue or cervix. Posteriorly, the Mesh Implant may be tacked to the pre-rectal tissue or cervix.

The VSD supports the vaginal tissues after surgery and facilitates abutment of the vaginal tissues against the Mesh Implant until tissue ingrowth occurs. Tissue ingrowth through the Mesh Implant occurs during the 3 to 4 weeks following surgery. Use of the GYNECARE PROSIMA System avoids the need for dissection outside the pelvic cavity, and avoids passage of suture and instruments through the obturator foramen and sacrospinous ligament, thus making surgery simpler to perform.

Hysterectomy

The surgeon's preference and the patient's needs determine if a concomitant hysterectomy is required. When a hysterectomy is performed, closure of the cul-de-sac peritoneum is recommended to avoid contact of the Mesh Implant with the bowel. A "T" incision closure should be avoided as this may increase the risk of mesh exposure. When vaginal hysterectomy is performed together with either or both anterior and posterior repair, the hysterectomy incision should first be closed transversely, and then the repair incisions should be made such that they do not connect with the previously closed hysterectomy incision. This is done to prevent the creation of a "T" incision.

Uterine Conservation

The GYNECARE PROSIMA System is suitable for use in situations when the surgeon or patient elects to conserve the uterus.

Vaginal Incisions

The vaginal incisions in the procedure using the GYNECARE PROSIMA System are the same as those used by the surgeon for routine vaginal repair surgery. Incisions should be made through the full depth of the vaginal wall to reduce the potential of mesh exposure.

Mesh Implant Placement

The Mesh Implants are held in place by the VSD until tissue ingrowth occurs. Therefore, it is unnecessary to fix the Mesh Implant straps into position. The apical portion of the Mesh Implant may be tacked onto the fascia in the midline at the vaginal apex using suture such as 2-0 MONOCRYL™ (Poliglecaprone 25) Suture or 2-0 Coated VICRYL™ (polyglactin 910) Suture. The vaginal epithelium should not be sutured onto the Mesh Implant.

Vaginal Preservation

Removing or excising too much vaginal epithelium should be avoided. Some tissue retraction may occur following surgery, and reduced vaginal capacity may be worsened if too much vaginal epithelium has been removed.

Three Levels of Vaginal Support

There are 3 levels of support to the vagina commonly known for vaginal repair. Use of the GYNECARE PROSIMA System in a procedure is intended to provide level I and II of this support as follows:

Level I - Suspension and Support (upper third of vagina)

The upper third of the vagina (including the vault following hysterectomy) and uterus are supported by 2 mechanisms. Firstly, direct support for the uterus and upper vagina is provided by the parametrium (cardinal and utero-sacral ligaments) and paracolpium fibers. These fibers act like suspensory ligaments and arise from the fascia of the piriformis muscle, sacrotuberous ligament, and lateral sacrum, and insert into the lateral upper third of the vagina and posterolateral aspect of the cervix. Secondly, indirect support for the uterus and upper vagina is provided by the levator plate, formed by the fusion of the right and left levator ani muscles between the rectum and coccyx. Uterine and vaginal vault prolapse occur as a consequence of failure of these direct and indirect supporting mechanisms. This is likely to involve weakness of the muscular pelvic floor and suspensory fibers of the parametrium and upper paracolpium. The aim of prolapse surgery at level I is to recreate direct and indirect supporting mechanisms. The GYNECARE PROSIMA System uses the Mesh Implant

straps to abut onto each obturator internus muscle and the overlying parietal fascia in the anterior vaginal repair, and Mesh Implant straps about the sacrospinous ligaments in the posterior vaginal repair. This provides direct support by suspension and indirect support by providing a broad area of Mesh Implant support for the upper vagina and uterus.

Level II - Lateral Attachment (mid third of vagina)

The mid-vagina is attached laterally and directly to the muscles on the pelvic sidewall by the arcus tendineus fasciae pelvis (ATFP). At this level, the anterior and posterior vaginal walls are stretched between right and left lateral attachments. At level II, prolapse repair aims to reattach the lateral mid-vagina onto the muscles of the pelvic sidewall. Central defects of the mid-vagina also require support at level II. Use of the GYNECARE PROSIMA System in a procedure recreates lateral attachment of the vagina onto the pelvic sidewall muscles and also provides central fascial reinforcement after tissue ingrowth.

Level III – Fusion (lower third of vagina)

NOTE: Dissection in this area is not required when using the GYNECARE PROSIMA System.

At level III, anteriorly the lower third of the vagina fuses with the perineal membrane and urethra. Posteriorly, the lower third of the vagina fuses with the perineal body and levator ani muscles. The tissues in this area are repaired without Mesh Implant, as the Mesh Implant is not intended to be used in the lower third of the vagina. The GYNECARE PROSIMA System does not address level III support defects, though they may be addressed by concomitant procedures such as perineorraphy.

SECTION 3: INSTRUCTIONS FOR USE

NOTE: Figures provided in the beginning of this document should be referenced when reading this section.

Surgical Preparation

Surgery performed with the GYNECARE PROSIMA System may be carried out under general or regional anesthesia according to the surgeon's, anesthesiologist's, and patient's preference.

The patient should be placed in the lithotomy position with the buttocks slightly overhanging the operating table and the hips flexed. At the surgeon's discretion the bladder may be drained. A catheter is required prior to balloon inflation and may be inserted at this point in the procedure.

GYNECARE PROSIMA System Use in Procedures Post-Hysterectomy

Anterior Vaginal Repair

When reinforcement of only the anterior vaginal wall is needed, only the GYNECARE PROSIMA Anterior Pelvic Floor Repair System should be used. This contains 1 Mesh Implant and a specially designed Anterior Inserter for use in an anterior vaginal repair. After the required vaginal incisions and dissections are made, tissue channels are created in the anterior compartment for placement of the Mesh Implant straps using the Anterior Inserter. **NOTE: The Anterior Inserter is not to be used for tissue dissection.**

Anterior Vaginal Dissection

The anterior vaginal epithelium is dissected off the bladder. Dissect the full thickness of the vaginal wall. This dissection should be facilitated by subepithelial hydrodissection. Superficial dissection of the vaginal wall or separating the vaginal wall into 2 layers should be avoided. Such dissection may result in a very thin vaginal wall and may also compromise vaginal wall blood supply, increasing the risk of mesh exposure. Laterally, continue the dissection towards the pelvic sidewall and to the ischial spine.

Anterior Channel Dissection and Mesh Implant Placement

For the purpose of this description, perform the dissection to create channels for the Mesh Implant straps first on the right side of the patient, then on the left. These channels are created in order to place the Mesh Implant such that the distal section of the straps lies flush against the pelvic side wall and parietal fascia of the obturator internus muscle. To place these straps, commence dissection by palpating and identifying the ischial spine on both sides. **NOTE: This dissection may alternatively be started with scissors using a "push-spread" technique, so that the tips of the scissors always remain anterior to the ischial spine.** Follow the initial dissection by gentle finger dissection to the ischial spine. Once contact is made with the ischial spine, sweep the index finger to create a space anterior and superior to the ischial spine. See Figure 8A. The direction of this dissection is perpendicular to the pelvic sidewall and creates a space approximately 2 cm in width and 3 cm in height. The anterior dissection does not involve dissection onto the sacrospinous ligaments. This dissection creates a channel anterior and superior to the ischial spine and superficial to the ATFP, the obturator internus muscle, and its parietal fascia. Repeat the same dissection on the left side.

Plication of the pre-vesical tissue is not required. However, if plication is performed, only the central portion of this tissue is plicated. This avoids making the dissected area too narrow. Place the Mesh Implant over the pre-vesical tissue with the strap pockets facing upwards. If tacking is to be performed it should be done at this point in the procedure by placing suture such as 2-0 MONOCRYL Suture or 2-0 Coated VICRYL Suture in the apex of the vagina and threading through the apical tab of the Mesh Implant. The tack may be tied down at this time or after the straps are placed. Tacking the distal groove of the Mesh Implant is optional and may be done with suture such as 2-0 MONOCRYL Suture or 2-0 Coated VICRYL Suture.

Using the Anterior Inserter, place the Mesh Implant straps into each right and left channel created by the dissection anterior and superior to the ischial spine (as described above). **NOTE: The curved ends of the Anterior Inserter are twisted in opposite directions, and there are arrows on each end indicating direction for placement.** With the arrow pointing towards the patient's right side, insert the tip of the Anterior Inserter into the Mesh Implant strap pocket (see Figure 8B) on the patient's right side. **NOTE: Counter-traction may help to keep the pocket loaded on the Anterior Inserter.** Keep the Anterior Inserter in a vertical position, such that the curved part of the instrument is against the posterior vaginal wall. Then direct the Anterior Inserter, with strap loaded, into the previously created tissue channel (see Figure 8C) until the handle comes in contact with the labia majora on the contra-lateral side. This is accomplished by positioning the handle portion of the Anterior Inserter in an upward-vertical direction such that the leading edge and pocket goes towards the ischial spine. Once positioned, cant the handle downward to the near-horizontal position, while maintaining the handle in contact with the contra-lateral thigh. **NOTE: Retracting the bladder with a standard surgical instrument may be helpful for initial placement in the channel. If desired, use an index finger in the channel to guide the initial placement of the Anterior Inserter against the labia majora on the contra-lateral side, prior to lowering the handle.** Pushing slightly upward ensures that the strap pockets are positioned appropriately and the apical section of the Mesh Implant will abut the vaginal apex. **NOTE: If resistance is felt during strap insertion, determine the cause before proceeding. Continuing to advance the inserter while under resistance may result in damage to the Mesh Implant, or over-insertion causing damage to critical tissue structures.**

To remove the Anterior Inserter, cant the handle back to the vertical position before withdrawing, leaving the strap in the channel. **NOTE: Insert first strap completely. NOTE: If the Anterior Inserter is pulled out before the Mesh Implant strap is delivered to the target, the strap will need to be removed, re-loaded, and re-inserted.** Repeat this on the opposite side of the patient by flipping the Anterior Inserter and inserting the end, with the arrow pointing to the patient's left side, into the other pocket. Figure 8D shows both straps placed. **NOTE: During placement of the second strap, take care to avoid movement of the Mesh Implant and confirm that the Mesh Implant is NOT twisted.**

Position the body of the Mesh Implant loosely over the underlying vaginal tissue. Folding or twisting of the body and straps should be avoided. The Mesh Implant body may require trimming depending on the vaginal dimensions or amount of lateral dissection. The vaginal epithelium may be trimmed, but excessive removal of vaginal epithelium

should be avoided. Close the epithelium over the Mesh Implant without using interlocking sutures (as described below, see Figure 8E). The final placement of the Mesh Implant in the anterior compartment is shown in Figure 8F.

NOTE: Ensure hemostasis is achieved before and during closure of the vaginal incisions.

Close the vaginal incisions without interlocking or figure-of-eight sutures. This is to avoid de-vascularizing the vaginal epithelium along the incision lines and to reduce mesh enson. Preferably the epithelium is closed in 2 layers to obtain a relatively thick suture line at the site of the vaginal incision. Close the deeper layer using a continuous subepithelial non-interlocking stitch with suture such as 2-0 MONOCRYL Suture or 2-0 MONOCRYL™ Plus Antibacterial (Poliglecaprone 25) Suture. Then close the epithelium by a non-interlocking continuous everting mattress stitch, using suture such as 2-0 Coated VICRYL Suture or 2-0 Coated VICRYL™ Plus (polyglactin 910) Antibacterial Suture. **NOTE: Place Mesh Implant in the upper 2/3 of vagina, taking care to trim Mesh Implant if beyond upper 2/3.** If not done already, cystoscopy is recommended to exclude urinary tract injury.

Alternatively, a single-layered closure of the vaginal wall can be performed. A continuous everting, non-interlocking mattress stitch or interrupted stitches of suture such as 2-0 Coated VICRYL Suture or 2-0 Coated VICRYL Plus Suture may be used.

Posterior Vaginal Repair

When reinforcement of only the posterior vaginal wall is needed, only use the GYNECARE PROSIMA Posterior Pelvic Floor Repair System. This contains 1 Mesh Implant and a specially designed Posterior Inserter that is used for the posterior vaginal repair. After making the required vaginal incisions and dissections, create tissue channels in the posterior compartment in which to place the Mesh Implant straps. **NOTE: The Posterior Inserter is not to be used for tissue dissection.**

Posterior Vaginal and Channel Dissection

Dissect the posterior vaginal epithelium off the pre-rectal tissue. As with the anterior vaginal wall, the full thickness of the posterior vaginal wall should be dissected. This dissection should be facilitated by subepithelial hydro-dissection. Continue the dissection laterally on each side to the levator ani muscles at the level of the ischial spine. Then continue dissection through each of the rectal pillars and onto, but not through, each sacrospinous ligament, creating channels in which the Mesh Implant straps will be placed. See Figure 9A.

Management of preexisting enterocele is optional, but if performed, may be carried out at this stage according to the surgeon's preferred technique.

If the peritoneal cavity is opened during either anterior or posterior dissection, it must be closed prior to mesh placement.

Posterior Mesh Implant Placement

Plication of the pre-rectal tissue is not required. However, if plication of the pre-rectal tissue is performed, only the central portion of the pre-rectal tissue is plicated. This avoids making the dissected area too narrow. Place the Mesh Implant over the pre-rectal tissue with the strap pockets facing upwards. If tacking is to be performed it should be done at this point in the procedure by placing suture such as 2-0 MONOCRYL Suture or 2-0 Coated VICRYL Suture in the apex of the vagina and threading through the apical tab of the Mesh Implant. The tack may be tied down at this time or after the straps are placed. Tacking the distal groove of the Mesh Implant is optional and may be done with suture such as 2-0 MONOCRYL Suture or 2-0 Coated VICRYL Suture.

Using the Posterior Inserter, place the Mesh Implant straps into each right and left channel created by the dissection towards each sacrospinous ligament (as described above). Grasp the Posterior Inserter using a straight needle holder/driver, as shown in Figure 9B. **NOTE: Place the tip of the needle holder/driver inside the straight grooved end of the Posterior Inserter.** Ensure that the connected Posterior Inserter is in-line with the handle of the needle holder/driver. Insert the tip of the Posterior Inserter into the strap pocket on the patient's right side (see Figure 9B). Then, direct the Posterior Inserter, with strap loaded, into the previously created tissue channel (see Figure 9C), keeping the position of the handle of the needle holder/driver upright. Proceed to insert the full length of the strap into the channel so that the base of the strap meets the superior limit of the fascial dissection. **NOTE: Insert first strap completely. If the inserter is pulled out before the strap is delivered to the target, the strap will need to be removed, re-loaded, and re-inserted.** **NOTE: Take care not to insert too deeply to avoid damage to critical tissue structures. NOTE: If resistance is felt during strap insertion, determine the cause before proceeding. Continuing to advance the inserter while under resistance may result in damage to the Mesh Implant or over-insertion, causing damage to critical tissue structures.** Withdraw the Posterior Inserter along the insertion path, leaving the strap in the channel. The straps abut, but do not penetrate, the sacrospinous ligaments. Do not place sutures into the sacrospinous ligaments. Repeat the procedure on the patient's left side with the second strap. Figure 9D shows both straps placed. **NOTE: During placement of the second strap, take care to avoid movement of the Mesh Implant and confirm that the Mesh Implant is NOT twisted.**

Position the body of the Mesh Implant loosely over the underlying vaginal fascia. Avoid folding or twisting the Mesh Implant body and straps. Depending on the vaginal dimensions, or the amount of lateral dissection, the Mesh Implant body may require trimming. The posterior vaginal wall epithelium may be trimmed but excessive removal of vaginal epithelium should be avoided. Close the posterior vaginal wall epithelium over the Mesh Implant without using interlocking sutures (as described below). The final placement of the Mesh Implant in the posterior compartment is shown in Figure 9E.

NOTE: Ensure hemostasis is achieved before and during closure of the vaginal incisions.

Close the vaginal incisions without using interlocking or figure-of-eight sutures. This is to avoid de-vascularizing the vaginal epithelium along the incision lines and to reduce mesh erosion. Preferably close the epithelium in 2 layers to obtain a relatively thick suture line at the site of the vaginal incision. Close the deeper layer using a continuous subepithelial non-interlocking stitch with suture such as 2-0 MONOCRYL Suture or 2-0 MONOCRYL Plus Antibacterial Suture. Then close the epithelium with a non-interlocking continuous everting mattress stitch, using suture such as 2-0 Coated VICRYL Suture or 2-0 Coated VICRYL Plus Suture. **NOTE: Place Mesh Implant in the upper 2/3 of vagina, taking care to trim Mesh Implant if beyond upper 2/3.** At the completion of surgery, a digital rectal examination is required to exclude rectal injury.

Alternatively, a single-layered closure of the vaginal wall can be performed. A continuous everting, non-interlocking mattress stitch or interrupted stitches of suture such as 2-0 Coated VICRYL Suture or 2-0 Coated VICRYL Plus Suture may be used.

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When both anterior and posterior vaginal wall reinforcement is needed, the GYNECARE PROSIMA Combined Pelvic Floor Repair System is used. This contains 2 identical Mesh Implants, one for the anterior vaginal repair and the second for the posterior vaginal repair. Use only the curved Anterior Inserter for an anterior repair and only the straight Posterior Inserter for a posterior repair. Perform the anterior and posterior vaginal repairs as described above. It is recommended that the anterior vaginal repair be performed first. The final placement of Mesh Implants in the anterior and posterior compartments is shown in Figure 10. At the completion of surgery, cystoscopy is recommended to exclude urinary tract injury. A digital rectal examination is required to exclude rectal injury.

GYNECARE PROSIMA System Use with Uterine Preservation (Hysteropexy)

If the prolapsed uterus is conserved, the apical tab of the Mesh Implant should be fixed to the cervix. Fixation of the Mesh Implant onto the cervix should occur at the level of the pubo-cervical ring when placed during the anterior or posterior vaginal repair.

When the uterus is conserved during an anterior vaginal repair the pubo-cervical ring is exposed during the anterior vaginal dissection. Place a 2-0 PROLENE Suture firmly into the anterior aspect of the pubo-cervical ring. This suture is also placed through the apical tab of the Mesh Implant. The PROLENE Suture at the tab is tied after the Mesh Implant straps are in place. This secures the Mesh Implant to the anterior surface of the cervix at the level of the pubo-cervical ring and ensures the Mesh Implant distends with the vagina as the VSD is properly positioned.

In the posterior repair, attach the Mesh Implant to the posterior cervix at, or above, the level of the pubo-cervical ring. The cul-de-sac may be opened during the attachment of the Mesh Implant to the cervix. Close the peritoneum of the cul-de-sac above this suture to prevent bowel adhering to the Mesh Implant. If the surgeon chooses not to open the cul-de-sac, the pubo-cervical ring is exposed during the posterior vaginal dissection. A 2-0 PROLENE Suture is placed firmly into the posterior aspect of the pubo-cervical ring. This suture is also placed through the apical tab of the Mesh Implant. The PROLENE Suture is tied after the Mesh Implant straps are in place. This secures the Mesh Implant to the posterior surface of the cervix at the level of the pubo-cervical ring.

When used for both anterior and posterior vaginal repairs, the Mesh Implants should be fixed to the anterior and posterior aspects of the cervix as described above (see Figure 11).

Mesh Implant Hygiene

During surgery, irrigate the vaginal wounds with saline. Keep the handling of the Mesh Implant to a minimum and practice good mesh hygiene.

Placement of the VSD and Balloon

At the completion of surgery, place an appropriately sized VSD with attached Balloon in the vagina and suture it in position to prevent dislodgement. The VSD has 3 potential sizes (small, medium and large) and can be customized by the surgeon to suit the vaginal length of the patient as follows.

VSD Fitting and Trimming

The VSD is provided in its largest size. Determine the appropriate size of the VSD for the patient by using the VSD itself to assess fit to the patient. This is done by placing the large size VSD in the vagina between the distended apex and the hymenal ring. To insert the VSD into the vagina, grasp at the widest point of the VSD and fold along the longitudinal axis with the Balloon facing upwards (see Figure 12). The widest point of the VSD is inserted first so that the suture holes are located just above the hymenal ring. **NOTE: Do not remove or damage the Balloon during VSD sizing.** Proper size is achieved when the VSD fits snugly in the upper 2/3 of the distended vagina with the distal end and suture eyelets 1 cm above the hymenal ring (see Figure 13).

If the large size fits then the VSD is not modified. If the medium size is required, then the uppermost section is removed by carefully trimming using only the tips of curved Mayo scissors to take small bites and to ensure a smooth cut edge. Care should be taken to minimize the amount of material remaining at the cut areas. **NOTE: It is important to fit the VSD very carefully. Once a VSD is cut it cannot be made larger and the cut sections cannot be reattached.** Move the Balloon out of the way during trimming (see Figure 14). **Care should be taken to avoid damaging the Balloon during trimming of the VSD.**

If the medium size fits then no further trimming is required. If the small size is required, then the remaining section is removed as above. Move the Balloon out of the way during trimming to avoid damaging it.

Once the VSD is properly sized and the Balloon is repositioned, the assembly can be inserted into the patient's vagina. **NOTE: To minimize risk of perforation of the Balloon, do not use any instruments to assist with insertion of the VSD or of the Balloon.** If the Balloon becomes damaged, remove the Balloon from the VSD and use gauze packing to fill the vaginal cavity.

After the assembly is properly positioned in the upper 2/3 of the patient's distended vagina, secure the VSD in place by placing a single throw of suture through each VSD suture eyelet and into the posterior vaginal wall epithelium, laterally and above the hymen on each side, as shown in Figure 15, at the 4 and 8 o'clock positions. The right and left sutures are then tied in turn, holding the VSD firmly in position within the vagina. **NOTE: Take care not to puncture the Balloon when suturing the VSD into position.** 2-0 Coated VICRYL Suture or equivalent absorbable suture is recommended for this application.

Balloon Inflation

After suturing the VSD into position, attach the provided 50-mL syringe by twisting to lock it onto the Balloon's valve. **NOTE: After placing the VSD, a catheter is required to avoid urinary retention.** Following inflation with a small volume of ambient air (see Figure 16), palpate the full length of the Balloon with a finger to ensure the Balloon has deployed and is seated to the full extent of the vagina. Once deployment is confirmed, remove finger and continue to inflate the Balloon fully until only a fingertip fits snugly at the introitus between the Balloon and the vaginal wall. Stabilization of the VSD is recommended as inflation occurs. The inflated Balloon serves to abut the Mesh Implant onto the vaginal wall. The volume of air required to sufficiently inflate the Balloon will vary from patient to patient. **NOTE: The maximum Balloon inflation volume must not exceed 90 mL.** Once adequately inflated, detach the syringe from the valve by twisting. The Balloon's inflation line must extend out of the vagina to be affixed to the patient's thigh. The cap must be secured to the Balloon's valve to ensure that the Balloon will maintain the intended volume of air (see Figure 7). **NOTE: Do not over-tighten the cap.** If needed, the Balloon can be adjusted later, using a standard syringe to increase or decrease the volume of air within the Balloon. At any time the Balloon may be palpated or visually inspected to ensure it has maintained sufficient inflation. **NOTE: As the patient moves, the Balloon will settle in the vaginal cavity and may appear to either increase or decrease in pressure. This is normal.**

NOTE: Do not detach Balloon from VSD prior to use.

NOTE: Do not inflate the Balloon prior to its insertion in the vagina.

NOTE: After Balloon inflation, if the VSD suture eyelets have moved more than 1 cm above the hymenal ring or if there is excessive tension on the eyelet sutures, then decrease pressure on the Balloon, and if needed, reposition or resize the VSD.

NOTE: If any holes are noticeable in the Balloon, or if a leak is detected, or if the Balloon fails to remain expanded after inflation, then DO NOT use the Balloon. It should be removed from the VSD and disposed of by proper means. Use standard gauze packing in place of the Balloon.

NOTE: If the Balloon's connector plug detaches from the VSD, it should be pushed back in place.

NOTE: Do not secure the Balloon inflation line in the vagina.

NOTE: In order to prevent damage, never apply extreme bending, tension, or twisting forces to the inflation line.

NOTE: Do not apply gauze packing in the presence of a Balloon.

Balloon Removal from VSD

Using a standard syringe, completely deflate and remove the Balloon 1 day after surgery, leaving the VSD in place. **NOTE: Do not leave the Balloon inside the vagina for more than 1 day.**

1) Remove the cap from the Balloon's valve.

2) Attach a standard 50-mL (or larger) syringe to the Balloon's valve and fully deflate the Balloon (see Figure 17). It is important to completely deflate the Balloon before attempting to remove it from the VSD. **NOTE: A fully deflated balloon will cause the syringe plunger to retract after removal of all air.**

3) Remove the syringe.

4) The Balloon can then be separated from the VSD and removed from the patient by gentle pulling in a caudal direction on the inflation line at a location near the Balloon's connector plug, while providing gentle counter-traction on the distal end of the VSD with a finger. See Figure 18.

NOTE: Do not retract the Balloon unless it is fully deflated and no resistance is felt. If resistance is felt, determine the cause before proceeding. Continuing to advance or retract the Balloon while under resistance may result in movement of the VSD and/or tissue trauma to the vaginal cavity. To ensure complete deflation has occurred, reconnect the syringe and remove all air before continuing with removal.

VSD Removal from Patient

Remove the VSD from the patient approximately 3 to 4 weeks following surgery after sufficient healing has occurred. By this time the absorbable sutures may have dissolved or lost sufficient tensile strength to allow easy removal of the VSD without any suture resistance. **NOTE: Cutting of both sutures may be required for removal. NOTE: Do not leave the VSD inside the vagina for more than 4 weeks.** Remove any remaining VSD attachment sutures. Manually remove the VSD from vaginal canal, as shown in Figure 19.

Perioperative Care

Patients can receive prophylactic antibiotics administered according to the surgeon's usual practice. Antibiotics may be continued postoperatively depending on the preference of the surgeon. Thromboembolic prophylaxis may be used.

The surgeon should explain the purpose of the VSD, which remains in the vagina for up to four weeks following surgery, is to support the vagina against the mesh during the healing period. The patient should be advised that the VSD will be removed during a post-operative check-up, approximately 4 weeks following surgery. The patient should be advised that postoperative vaginal discharge may be experienced and that the VSD may move down slightly. If the patient feels that the VSD has moved down, she may gently push it up to a more comfortable position. However, if the VSD is causing significant discomfort, the patient should be informed to contact their doctor.

Following discharge from the hospital, the patient should be instructed to avoid strenuous activity for a period of 3 to 4 weeks. By this time the pelvic tissues will have incorporated into the Mesh Implant, and the patient can then resume activities of normal daily living. The patient should be advised to avoid sexual intercourse for at least 6 weeks following surgery. Pelvic floor exercises may be recommended any time after surgery.

PERFORMANCE

Animal studies show that implantation of GYNECARE GYNEMESH PS elicits a minimal to slight inflammatory reaction, which is transient and is followed by the deposition of a thin fibrous layer of tissue which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The mesh remains soft and pliable, and normal wound healing is not noticeably impaired. The material is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes.

CONTRAINDICATIONS

- When GYNECARE GYNEMESH PS is used in infants, children, pregnant women, or women planning future pregnancies, the surgeon should be aware that this product will not stretch significantly as the patient grows.
- The GYNECARE PROSIMA System should not be used in the presence of pregnancy or purulent infections or cancers of the vagina, cervix, or uterus.

WARNINGS AND PRECAUTIONS

- Users should be familiar with surgical procedures and techniques involving pelvic floor repair and nonabsorbable meshes before employing the GYNECARE PROSIMA Systems.
- Use of the GYNECARE PROSIMA System has not been fully evaluated in patients with Stage IV pelvic organ prolapse. Therefore its use in these patients is not recommended.
- Acceptable surgical practice should be followed for the GYNECARE PROSIMA System as well as for the management of infected or contaminated wounds.
- Do not use the GYNECARE PROSIMA System if you think the surgical site may be infected or contaminated. If the Mesh Implant or VSD-Balloon Assembly is used in contaminated areas it must only be with the understanding that subsequent infection may require its removal.
- Postoperatively the patient should be advised to refrain from heavy lifting and/or exercise (e.g. cycling, jogging) for 3 to 4 weeks and to refrain from sexual intercourse for 6 weeks or until the physician determines it is suitable for the patient to return to her normal activities.
- Do not leave the VSD inside the vagina for longer than 4 weeks.
- Do not leave the Balloon inside the vagina for longer than 1 day.
- The GYNECARE PROSIMA System components are not intended to be used with devices other than those mentioned in this package insert.
- Avoid placing excessive tension on the Mesh Implant during handling.
- Use the GYNECARE PROSIMA Systems with care, and with attention to patient anatomy, to avoid damage to vessels, nerves, bladder, bowel, and vaginal wall perforation. Correct use of the GYNECARE PROSIMA System components will minimize risks.
- Inflate the Balloon only with ambient air.
- Palpation will confirm that the Balloon does not contain any air leaks after inflation. Complete loss of inflation may limit the Balloon's effectiveness.
- The Balloon wall is thin in order to achieve desired properties. Punctures, cuts, nicks, crushing, or overstretching can lead to a loss of inflation. The Balloon may be easily penetrated by needle or scalpel or ruptured by manipulation with a blunt instrument. Care must be exercised during handling to prevent such events. A damaged Balloon must not be used. Remove and pack with gauze.
- The Balloon inflation maximum is 90 mL. Do not over-inflate the Balloon. Excessive inflation of the Balloon may cause patient discomfort, tissue necrosis, disruption of vaginal wound postoperatively, or inability to void.
- Do not use GYNECARE PROSIMA Systems on patients who are on anti-coagulant therapy.

- Bleeding may occur postoperatively. Observe for any symptoms or signs before releasing the patient from the hospital.
- The patient should be instructed to contact the surgeon immediately if unusual pain, bleeding, or other problems occur.
- Although bladder injury is unlikely to occur with this technique, cystoscopy is recommended to be performed.
- Although rectal injury is unlikely to occur with this technique, a digital exam is required to be performed.
- Do not affix the GYNECARE GYNEMESH PS Mesh Implant with any staples, clips, or clamps as mechanical damage to the mesh may occur.
- The Mesh Implant should not be present in the lower 1/3 of vagina. If needed, trim the Mesh Implant to the junction of the lower and middle 1/3 of vaginal wall.
- Prophylactic antibiotics can be administered according to the surgeon's usual practice.

ADVERSE REACTIONS

- Potential adverse reactions are those typically associated with surgically implantable materials, including infection potentiation, inflammation, adhesion formation, fistula formation, erosion, extrusion and scarring that result in implant contraction.
- Potential adverse reactions are those typically associated with pelvic organ prolapse repair procedures, including pain with intercourse and pelvic pain. These may be self-resolving over time.
- Punctures or lacerations or injury to vessels, nerves, bladder, urethra, or bowel may occur during dissection or mesh placement and may require surgical repair.
- Dissection for pelvic floor repair procedures has the potential to impair normal voiding for a variable length of time.

STERILITY

The GYNECARE PROSIMA Systems are sterilized by ethylene oxide. DO NOT RESTERILIZE any portion of the GYNECARE PROSIMA System. DO NOT REUSE any portion of the GYNECARE PROSIMA System. Reuse of this device (or portions of this device) may create a risk of product degradation and cross-contamination, which may lead to infection or transmission of bloodborne pathogens to patients and users. Do not use if package is opened or damaged. Discard all opened, unused GYNECARE PROSIMA System components.

DISPOSAL

Dispose of the GYNECARE PROSIMA System components and packaging according to your facility's policies and procedures concerning biohazardous materials and waste.

STORAGE

Recommended storage conditions: controlled room temperature and relative humidity (approximately 25°C, 60% RH), away from moisture and direct heat. Do not use after expiry date.

Symbols Used on Labeling

 0086 CE mark and identification number of Notified Body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.	 Manufacturer
 Batch number	 Do not reuse/resterilize
 Use by — year and month	 See instructions for use  Method of Sterilization — Ethylene Oxide

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Johnson & Johnson Medical Limited
CH Johnson & Johnson AG