

EXHIBIT Z

HISTION**Accession Log****DOC010**

Study Number	H16-008
---------------------	---------

Accession #	Animal #	Group #	Description
16-008- 1	157162		Paraffin
16-008- 2	157164		Paraffin
16-008- 3	157166		Paraffin
16-008- 4	157170		Paraffin
16-008- 5	157172		Paraffin
16-008- 6	157174		Paraffin
16-008- 7	157178		Paraffin
16-008- 8	157180		Paraffin
16-008- 9	157182		Paraffin
16-008- 10	157186		Paraffin
16-008- 11	157188		Paraffin
16-008- 12	157190		Paraffin
16-008- 13	157194		Paraffin
16-008- 14	157198		Paraffin
16-008- 15	157232		Paraffin
16-008- 16	157234		Paraffin
16-008- 17	157236		Paraffin
16-008- 18	157880		Paraffin
16-008- 19	157884		Paraffin
16-008- 20	157895		Paraffin
16-008- 21	157899		Paraffin
16-008- 22	157903		Paraffin
16-008- 23	159703		Paraffin
16-008- 24	159704		Technovit
16-008- 25	157163		Technovit
16-008- 26	157165		Technovit
16-008- 27	157167		Technovit
16-008- 28	157169		Technovit
16-008- 29	157171		Technovit
16-008- 30	157173		Technovit
16-008- 31	157175		Technovit
16-008- 32	157177		Technovit
16-008- 33	157179		Technovit
16-008- 34	157181		Technovit
16-008- 35	157183		Technovit
16-008- 36	157185		Technovit
16-008- 37	157187		Technovit
16-008- 38	157189		Technovit
16-008- 39	157191		Technovit
16-008- 40	157193		Technovit
16-008- 41	157195		Technovit
16-008- 42	157199		Technovit
16-008- 43	157201		Technovit
16-008- 44	157233		Technovit
16-008- 45	157235		Technovit
16-008- 46	157237		Technovit
16-008- 47	Rabbit Skin		Paraffin
16-008- 48	Rabbit Skin		Technovit

Checked By:



Date: 3/10/16

Confidential

H16-008_study forms paraffin_08Mar16_TAC v3.2

Histology Grossing Log

DOC005

Study Number		H16-008	
Specimen Type		Mesh	
Grossing SOP/SSP		N/A	
Description of grossing method		Drawing	
<p>Specimens cut along lines marked on attached photographs; allowing for trimming in to plane of section (see drawing) and to prevent unravelling.</p>			
After grossing, specimens transferred to:			
ID	A	B	C
Dissector	JDA		
Date	3/10/16		
Enter the name of the dissector and the date for each grossing day. When tissue grossing complete, enter the ID letter in the Done box. If there are no observations, enter N/A.			
Specimen #	Done?	Observations/comments/deviations	Block IDs
16-008-1	A	N/A	
16-008-2	A	N/A	
16-008-3	A	N/A	
16-008-4	A	N/A	
16-008-5	A	N/A	
16-008-6	A	N/A	
16-008-7	A	N/A	
16-008-8	A	N/A	
16-008-9	A	N/A	
16-008-10	A	N/A	
16-008-11	A	N/A	
16-008-12	A	N/A	
16-008-13	A	N/A	
16-008-14	A	N/A	
16-008-15	A	N/A	
16-008-16	A	Lost specimen when trimming - processed other half	Ret *
16-008-17	A	N/A	
16-008-18	A	N/A	
16-008-19	A	N/A	
16-008-20	A	N/A	
16-008-21	A	N/A	
16-008-22	A	N/A	
16-008-23	A	N/A	
16-008-24	A	N/A	
16-008-25	A	N/A	
16-008-26	A	N/A	
16-008-27	A	N/A	
16-008-28	A	N/A	
16-008-29	A	N/A	
16-008-30	A	N/A	
16-008-31	A	N/A	
16-008-32	A	N/A	
16-008-33	A	N/A	
16-008-34	A	N/A	

* other piece (lost piece) found and returned to original package unprocessed JDA 3/13/16

DOC005

[illegible]

Processing Log							DOC033
Study Number	H16-008						
Specimen Type	Mesh						
Embedded In	Paraffin						
Specimen Numbers							
16-008--1	16-008--17	-	-	-	-	-	
16-008--2	16-008--18	-	-	-	-	-	
16-008--3	16-008--19	-	-	-	-	-	
16-008--4	16-008--20	-	-	-	-	-	
16-008--5	16-008--21	-	-	-	-	-	
16-008--6	16-008--22	-	-	-	-	-	
16-008--7	16-008--23	-	-	-	-	-	
16-008--8	16-008-47	-	-	-	-	-	
16-008--9	-	-	-	-	-	-	
16-008--10	-	-	-	-	-	-	
16-008--11	-	-	-	-	-	-	
16-008--12	-	-	-	-	-	-	
16-008--13	-	-	-	-	-	-	
16-008--14	-	-	-	-	-	-	
16-008--15	-	-	-	-	-	-	
16-008--16	-	-	-	-	-	-	
If processing using an automated tissue processor enter the run information in the table below.							
Processing			Paraffin Embedding				
EQP#	H15T035	Initial	EQP	H15T173			
SOP#	Protocol		Reagent	15-118			
Date Started	3/10/16		Date	3/11/16			
Date Finished	3/11/16		Initial				
If processing manually enter the run information in the table below.							
Duration = minimum required incubation time. Enter time when incubation steps are less than one day							
Manual Processing Steps							
Step #	Reagent Name	Duration	Reagent #	Date	Time (optional)	Initial	
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							

Processing Log							DOC033																																																																																																																																																			
Study Number	H16-008																																																																																																																																																									
Specimen Type	Mesh																																																																																																																																																									
Embedded In	Technovit																																																																																																																																																									
Specimen Numbers																																																																																																																																																										
16-008--24	16-008--40	-	-	-	-	-																																																																																																																																																				
16-008--25	16-008--41	-	-	-	-	-																																																																																																																																																				
16-008--26	16-008--42	-	-	-	-	-																																																																																																																																																				
16-008--27	16-008--43	-	-	-	-	-																																																																																																																																																				
16-008--28	16-008--44	-	-	-	-	-																																																																																																																																																				
16-008--29	16-008--45	-	-	-	-	-																																																																																																																																																				
16-008--30	16-008--46	-	-	-	-	-																																																																																																																																																				
16-008--31	-	-	-	-	-	-																																																																																																																																																				
16-008--32	-	-	-	-	-	-																																																																																																																																																				
16-008--33	-	-	-	-	-	-																																																																																																																																																				
16-008--34	-	-	-	-	-	-																																																																																																																																																				
16-008--35	-	-	-	-	-	-																																																																																																																																																				
16-008--36	-	-	-	-	-	-																																																																																																																																																				
16-008--37	-	-	-	-	-	-																																																																																																																																																				
16-008--38	-	-	-	-	-	-																																																																																																																																																				
16-008--39	-	-	-	-	-	-																																																																																																																																																				
<p><i>If processing using an automated tissue processor enter the run information in the table below.</i></p> <div style="display: flex; justify-content: space-around;"> <table border="1" style="width: 45%;"> <thead> <tr> <th colspan="2">Processing</th> <th>Initial</th> </tr> </thead> <tbody> <tr> <td>EQP#</td> <td>HISTO35</td> <td>pet</td> </tr> <tr> <td>SOP#</td> <td>Protocol</td> <td>pet</td> </tr> <tr> <td>Date Started</td> <td>3/10/16</td> <td>pet</td> </tr> <tr> <td>Date Finished</td> <td>3/10/16</td> <td>pet</td> </tr> </tbody> </table> <table border="1" style="width: 45%;"> <thead> <tr> <th colspan="2">Paraffin Embedding</th> </tr> </thead> <tbody> <tr> <td>EQP</td> <td rowspan="4" style="text-align: center; vertical-align: middle;">Z</td> </tr> <tr> <td>Reagent</td> </tr> <tr> <td>Date</td> </tr> <tr> <td>Initial</td> </tr> </tbody> </table> </div> <p><i>If processing manually enter the run information in the table below.</i></p> <p><i>Duration = minimum required incubation time. Enter time when incubation steps are less than one day</i></p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="7" style="text-align: center;">Manual Processing Steps</th> </tr> <tr> <th>Step #</th> <th>Reagent Name</th> <th>Duration</th> <th>Reagent #</th> <th>Date</th> <th>Time (optional)</th> <th>Initial</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Technovit 7200</td> <td>3hrs</td> <td>16-039</td> <td>3/10/16</td> <td>10pm</td> <td>pet</td> </tr> <tr> <td>2</td> <td>Technovit 7200</td> <td>3hrs</td> <td>16-039</td> <td>3/11/16</td> <td>1am</td> <td>pet</td> </tr> <tr> <td>3</td> <td>Technovit 7200</td> <td>3hrs</td> <td>16-039</td> <td>3/11/16</td> <td>4am</td> <td>pet</td> </tr> <tr> <td>4</td> <td>Polymerization</td> <td>4hrs/10hrs</td> <td>—</td> <td>3/11/16</td> <td>10am</td> <td>pet</td> </tr> <tr><td>5</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>6</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>7</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>8</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>9</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>10</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>11</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>12</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>13</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>14</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>15</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr><td>16</td><td></td><td></td><td></td><td></td><td></td><td></td></tr> </tbody> </table>							Processing		Initial	EQP#	HISTO35	pet	SOP#	Protocol	pet	Date Started	3/10/16	pet	Date Finished	3/10/16	pet	Paraffin Embedding		EQP	Z	Reagent	Date	Initial	Manual Processing Steps							Step #	Reagent Name	Duration	Reagent #	Date	Time (optional)	Initial	1	Technovit 7200	3hrs	16-039	3/10/16	10pm	pet	2	Technovit 7200	3hrs	16-039	3/11/16	1am	pet	3	Technovit 7200	3hrs	16-039	3/11/16	4am	pet	4	Polymerization	4hrs/10hrs	—	3/11/16	10am	pet	5							6							7							8							9							10							11							12							13							14							15							16						
Processing		Initial																																																																																																																																																								
EQP#	HISTO35	pet																																																																																																																																																								
SOP#	Protocol	pet																																																																																																																																																								
Date Started	3/10/16	pet																																																																																																																																																								
Date Finished	3/10/16	pet																																																																																																																																																								
Paraffin Embedding																																																																																																																																																										
EQP	Z																																																																																																																																																									
Reagent																																																																																																																																																										
Date																																																																																																																																																										
Initial																																																																																																																																																										
Manual Processing Steps																																																																																																																																																										
Step #	Reagent Name	Duration	Reagent #	Date	Time (optional)	Initial																																																																																																																																																				
1	Technovit 7200	3hrs	16-039	3/10/16	10pm	pet																																																																																																																																																				
2	Technovit 7200	3hrs	16-039	3/11/16	1am	pet																																																																																																																																																				
3	Technovit 7200	3hrs	16-039	3/11/16	4am	pet																																																																																																																																																				
4	Polymerization	4hrs/10hrs	—	3/11/16	10am	pet																																																																																																																																																				
5																																																																																																																																																										
6																																																																																																																																																										
7																																																																																																																																																										
8																																																																																																																																																										
9																																																																																																																																																										
10																																																																																																																																																										
11																																																																																																																																																										
12																																																																																																																																																										
13																																																																																																																																																										
14																																																																																																																																																										
15																																																																																																																																																										
16																																																																																																																																																										

HISTION

Staining Log

DOC008

Date	3/14/16
Initial	JS

QC date	3/14/16
Initial	JS

Slide#	QC	Slide#	QC	Slide#	QC	Slide#	QC
16-008-24	S	16-008-36	S				
16-008-25	S	16-008-37	S				
16-008-26	S	16-008-38	S				
16-008-27	S	16-008-39	S				
16-008-28	S	16-008-40	S				
16-008-29	S	16-008-41	S				
16-008-30	S	16-008-42	S				
16-008-31	S	16-008-43	S				
16-008-32	S	16-008-45	S				
16-008-33	S	16-008-46	S				
16-008-34	S	16-008-48	S				
16-008-35	S	16-008- written in error JS					

QC Codes

✓	Pass	K	Knife marks	T	Tissue missing	S	Stain Quality
x	Fail	D	Deeper	F	Fiber	μ	Thickness
W	Wrinkle/fold	C	Chatter	P	Pickup		
R	Recut	H	Trimming Holes	B	Bubbles		

Technovit plastic surrounding specimens stained too dark E eosin JS 3/14/16

Stain	H+E	SOP#	Protocol	Revision #	N/A
Manual	✓	Autostainer	EQP#	N/A	

If staining manually, please fill in the table below. If using the autostainer attach the process record to this sheet.

Record time in minutes. To indicate seconds use "s", hours use "hr".

	Reagent Name	Reagent #	Time	Comment	Done?
1	100% Alcohol	16-047	30s		✓
2	Distilled water	15-161 JS 3/14/16	1		✓
3	Harris Haematoxylin	15-161	10 mins		✓
4	Distilled water	R16-043 JS 3/14/16			✓
5	1% Acid alcohol	R16-043	30s		✓
6	Tapwater	—	1		✓
7	0.25% NH ₄ OH	R16-056	1		✓
8	Tapwater	—	1		✓
9	1% Aqueous Eosin	R15-300	1		✓
10	Distilled water	—	30s		✓
11	Blot + Air dry	—	—		✓
12	Mount in Technovit 7200 VLC	16-040	—		✓
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					

HISTION

Staining Log

DOC008

Date	3/15/16
Initial	fat

QC date	3/15/16
Initial	fat

Slide#	QC	Slide#	QC	Slide#	QC	Slide#	QC
16-008-24	✓	16-008-36	✓				
16-008-25	✓	16-008-37	✓				
16-008-26	✓	16-008-38	✓				
16-008-27	✓	16-008-39	✓				
16-008-28	✓	16-008-40	✓				
16-008-29	✓	16-008-41	✓				
16-008-30	✓	16-008-42	✓				
16-008-31	✓	16-008-43	✓				
16-008-32	✓	16-008-44	✓				
16-008-33	✓	16-008-45	✓				
16-008-34	✓	16-008-46	✓				
16-008-35	✓						

QC Codes

✓	Pass	K	Knife marks	T	Tissue missing	S	Stain Quality
x	Fail	D	Deeper	F	Fiber	μ	Thickness
W	Wrinkle/fold	C	Chatter	P	Pickup		
R	Recut	H	Trimming Holes	B	Bubbles		

Stain	H+E	SOP#	Protocol.	Revision #	NA
Manual	✓	Autostainer	EQP#	N/A	

If staining manually, please fill in the table below. If using the autostainer attach the process record to this sheet.

Record time in minutes. To indicate seconds use "s", hours use "hr".

	Reagent Name	Reagent #	Time	Comment	Done?
1	Distilled water	—	1 min		✓
2	Harris Haematoxylin	15-161	4 hrs.	Checked microscopically @ 15 min intervals	✓
3	Distilled water	—	1		✓
4	1% Acid alcohol	R16-043	30s		✓
5	Tapwater	—	1		✓
6	0.25% NH ₄ OH	R16-056	1		✓
7	Tapwater	—	1		✓
8	1% Eosin	R15-300	1		✓
9	Distilled water	—	30s		✓
10	Blot + air dry	—	—		✓
11	Technovit 7200 VLC	16-040			✓
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					

HISTION

Sectioning & Staining Assignment Log

DOC007

Study Number	H16-008								
Specimen Type	Mesh								
Embedded In	Paraffin								
Use the letters in the following table to identify the user for each day of sectioning									
ID	A	B	C	D	E	F	G	H	I
Date	3/11/16								
Equipment	HST099								
Initial	JS								
Thickness (µm)	5								
ID	J	K	L	M	N	O	P	Q	R
Date									
Equipment									
Initial									
Thickness (µm)									
ID	S	T	U	V	W	X	Y	Z	AA
Date									
Equipment									
Initial									
Thickness (µm)									
ID	BB	CC	DD	EE	FF	GG	HH	II	JJ
Date									
Equipment									
Initial									
Thickness (µm)									

Staining Procedure Table

Record the date of each staining event in the table below before completing the staining log.

Stain	H+E				
Initial	JS				
Date	3/12/16				
Stain					
Initial					
Date					
Stain					
Initial					
Date					
Stain					
Initial					
Date					
Stain					
Initial					
Date					
Stain					
Initial					
Date					

DOC007

[illegible]

HISTION

Staining Log

DOC008

Date	3/12/16
Initial	<i>[Signature]</i>

QC date	3/12/16
Initial	<i>[Signature]</i>

Slide#	QC	Slide#	QC	Slide#	QC	Slide#	QC
16-008-2	A X	16-008-6	A X				
16-008-8	A X	16-008-7	A X				
16-008-1	A X	16-008-9	A X				
16-008-3	A X	16-008-10	A X				
16-008-20	A X	16-008-11	A X				
16-008-21	A X	16-008-12	A X				
16-008-18	A X	16-008-13	A X				
16-008-22	A X	16-008-14	A X				
16-008-41	A X						
16-008-23	A X						
16-008-4	A X						
16-008-5	A X						

QC Codes All slides failed QC - Due to overstaining Eosin *[Signature]* 3/12/16

✓	Pass	K	Knife marks	T	Tissue missing	S	Stain Quality
x	Fail	D	Deeper	F	Fiber	μ	Thickness
W	Wrinkle/fold	C	Chatter	P	Pickup		
R	Recut	H	Trimming Holes	B	Bubbles		

Too Dark *[Signature]* 3/12/16

Stain	H+E	SOP#	Per Protocol	Revision #	
Manual	Autostainer	✓	EQP#	HIST100	

If staining manually, please fill in the table below. If using the autostainer attach the process record to this sheet.

Record time in minutes. To indicate seconds use "s", hours use "hr".

Timer

	Reagent Name	Reagent #	Time	Comment	Done?
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					

----- Process Report -----

Exported Date: 2016-03-12 10:50:27

Start:2016-03-12 10:03:50

End :2016-03-12 10:50:03

Program Name: he expo no 70

Specimen No : 1

Step	Station	Solution	Time	Result	Delay	Mix
1	S*	Start Station	--:--:--	00:00:15		
2	D*	Drying Station	00:10:00	00:10:00	**	OFF
3	3	Xylene	00:03:00	00:03:00	**	ON
4	5	Xylene	00:02:00	00:02:00	**	ON
5	6	Xylene	00:02:00	00:02:00	**	ON
6	17	Alcohol 100%	00:01:00	00:01:00	**	ON
7	15	Alcohol 100%	00:01:00	00:01:00	**	ON
8	31	Alcohol 95%	00:01:00	00:01:00	**	ON
9	W*	Wash Station	00:01:00	00:01:00	**	ON
10	22	Hematoxylin(Harris)	00:10:00	00:10:00	==	ON
11	W*	Wash Station	00:01:00	00:01:00	==	ON
12	23	Acid Alcohol	00:00:30	00:00:30	==	ON
13	W*	Wash Station	00:02:00	00:02:00	==	ON
14	24	Ammonia Water	00:01:00	00:01:17*	==	ON
15	W*	Wash Station	00:01:00	00:01:00	**	ON
16	26	Eosin	00:02:00	00:02:00	==	ON
17	16	Alcohol 100%	00:01:00	00:01:00	==	ON
18	18	Alcohol 100%	00:01:00	00:01:00	==	ON
19	19	Alcohol 100%	00:01:00	00:01:00	==	ON
20	6	Xylene	00:01:00	00:01:00	==	ON
21	5	Xylene	00:01:00	00:01:00	==	ON
22	E*	End Station	--:--:--	00:00:23		

----- Executed Diagnose -----

Step	Date	Description
13	2016-03-12 10:39:25	A deviation error occurred with the Z axis.
14	2016-03-12 10:40:18	Continue the process.

* Lid left on container- causing Z axis crash. Removed lid, placed slides in correct reagent and followed on-screen instructions. Slides spent 17 s extra in ammonia water. This had no effect on staining quality. JPT 3/12/16

DOC008

QC date	3/72/16
Initial	<i>[Signature]</i>

[illegible]

✓	Pass	K	Knife marks	T	Tissue missing	S	Stain Quality
✗	Fail	D	Deeper	F	Fiber	μ	Thickness
W	Wrinkle/fold	C	Chatter	P	Pickup		
R	Recut	H	Trimming Holes	B	Bubbles		

700 Dec.

Stain			SOP#	Per Protocol	Revision #
Manual		Autostainer	✓	EQP#	HIST 100

Timer		

Record time in minutes. To indicate seconds use "s", hours use "hr".

Timer		

----- Process Report -----

Exported Date: 2016-03-12 11:09:33

Start:2016-03-12 10:03:50

End :2016-03-12 11:05:14

Program Name: he expo no 70

Specimen No : 2

Step	Station	Solution	Time	Result	Delay	Mix
1	S*	Start Station	--:--:--	00:00:31		
2	D*	Drying Station	00:10:00	00:13:11	**	OFF
3	3	Xylene	00:03:00	00:03:00	**	ON
4	5	Xylene	00:02:00	00:02:00	**	ON
5	6	Xylene	00:02:00	00:02:05	**	ON
6	17	Alcohol 100%	00:01:00	00:01:00	**	ON
7	15	Alcohol 100%	00:01:00	00:01:00	**	ON
8	31	Alcohol 95%	00:01:00	00:13:44	**	ON
9	W*	Wash Station	00:01:00	00:01:00	**	ON
10	22	Hematoxylin(Harris)	00:10:00	00:10:00	==	ON
11	W*	Wash Station	00:01:00	00:01:00	==	ON
12	23	Acid Alcohol	00:00:30	00:00:30	==	ON
13	W*	Wash Station	00:02:00	00:02:00	==	ON
14	24	Ammonia Water	00:01:00	00:01:00	==	ON
15	W*	Wash Station	00:01:00	00:01:00	**	ON
16	26	Eosin	00:02:00	00:02:00	==	ON
17	16	Alcohol 100%	00:01:00	00:01:00	==	ON
18	18	Alcohol 100%	00:01:00	00:01:00	==	ON
19	19	Alcohol 100%	00:01:00	00:01:00	==	ON
20	6	Xylene	00:01:00	00:01:00	==	ON
21	5	Xylene	00:01:00	00:01:00	==	ON
22	E*	End Station	--:--:--	00:04:18		

----- Executed Diagnose -----

Step	Date	Description
8	2016-03-12 10:39:25	A deviation error occurred with the Z axis.
8	2016-03-12 10:40:18	Continue the process.

Z axis crash did not affect this set of slides *JS* 3/12/16

----- Stain Program List -----

Exported Date: 2016-03-12 09:43:57

Created Date : 2015-08-10 16:48:23

Revised Date : 2016-03-12 09:42:12

Program Name:he expo no 70

Step	Station	Solution	Time	Delay	Mix
1	S*	Start Station	--:--:--		
2	D*	Drying Station	00:10:00	**	OFF
3	3	Xylene	00:03:00	**	ON
4	5	Xylene	00:02:00	**	ON
5	6	Xylene	00:02:00	**	ON
6	17	Alcohol 100%	00:01:00	**	ON
7	15	Alcohol 100%	00:01:00	**	ON
8	31	Alcohol 95%	00:01:00	**	ON
9	W*	Wash Station	00:01:00	**	ON
10	22	Hematoxylin(Harris)	00:10:00	==	ON
11	W*	Wash Station	00:01:00	==	ON
12	23	Acid Alcohol	00:00:30	==	ON
13	W*	Wash Station	00:02:00	==	ON
14	24	Ammonia Water	00:01:00	==	ON
15	W*	Wash Station	00:01:00	**	ON
16	26	Eosin	00:02:00	==	ON
17	16	Alcohol 100%	00:01:00	==	ON
18	18	Alcohol 100%	00:01:00	==	ON
19	19	Alcohol 100%	00:01:00	==	ON
20	6	Xylene	00:01:00	==	ON
21	5	Xylene	00:01:00	==	ON
22	E*	End Station	--:--:--		

Total 22 Steps

HISTION

Paraffin sections - Hand stained

Staining Log

DOC008

Date	3/12/16
Initial	<i>[Signature]</i>

QC date	10/12/16
Initial	<i>[Signature]</i>

Slide#	QC	Slide#	QC	Slide#	QC	Slide#	QC
16-008-1	✓	16-008-13	✓	16-008-23	✓		
16-008-2	✓	16-008-14	✓				
16-008-3	✓	16-008-15	✓				
16-008-4	✓	16-008-16	✓				
16-008-5	✓	16-008-17	✓				
16-008-6	✓	16-008-18	✓				
16-008-7	✓	16-008-19	✓				
16-008-8	✓	16-008-47	✓				
16-008-9	✓	RACK 2.					
16-008-10	✓	16-008-20	✓				
16-008-11	✓	16-008-21	✓				
16-008-12	✓	16-008-22	✓				

QC Codes

✓	Pass	K	Knife marks	T	Tissue missing	S	Stain Quality
*	Fail	D	Deeper	F	Fiber	μ	Thickness
W	Wrinkle/fold	C	Chatter	P	Pickup		
R	Recut	H	Trimming Holes	B	Bubbles		

Stain	H+E	SOP#	Per Protocol.	Revision #
Manual	✓	Autostainer	EQP#	N/A

If staining manually, please fill in the table below. If using the autostainer attach the process record to this sheet.

Record time in minutes. To indicate seconds use "s", hours use "hr".

	Reagent Name	Reagent #	Time	Comment	Done?
1	Xylenes	16-005	3		✓
2	Xylenes	16-005	32		✓
3	Xylenes	16-005	2		✓
4	100% Alcohol	16-047	1		✓
5	100% Alcohol	16-047	1		✓
6	95% Alcohol	R16-042	1		✓
7	Distilled water.	-	1		✓
8	Morris Haematoxylin	15-161	10		✓
9	Tapwater	-	1		✓
10	Acid Alcohol	R16-043	30s		✓
11	Tapwater	-	2		✓
12	0.25% H ₂ O ₂	R16-056	1		✓
13	Tapwater	-	1		✓
14	1% Eosin	R15-300	2		✓
15	Tapwater	-		20s - Microscope differ	✓
16	70% Alcohol	R16-040	Rinse		✓
17	100% Alcohol	16-047	1		✓
18	100% Alcohol	16-047	1		✓
19	100% Alcohol	16-047	1		✓
20	Xylenes	16-005	2		✓
21	Xylenes	16-005	2		✓
22	Curemant II				✓
23					
24					
25					
26					
27					